

**"SEE & TREAT
IN JUST 3 mm!"**

GYNEWorld®
Women's Health

DIAGNOSTIC HYSTEROSCOPY



**Adaptable to optics:
0° / 12° / 30°**



**Up to 7 FR
working channel**

3 mm

7 FR CHANNEL

Distributed by:



Gynko®

**EXCLUSIVE WORLDWIDE
PATENTED PRODUCT**

info@gyneworld.com



ADVANTAGES

- ✓ TOTAL STERILITY IN OPTICS
- ✓ 3 mm IN DIAMETER FOR DIAGNOSTIC AND OPERATIVE PROCEDURES
- ✓ MULTIPLE PROCEDURES WITH A SINGLE OPTICAL DEVICE
- ✓ IT INCLUDES PROTECTIVE CAMERA CASE

GFRGYN300

Gynko® sheath with 7 FR working channel for 2.9 mm I=300 mm rigid optics, 0° with Storz connection.

GFRGYN330

Gynko® sheath with 7 FR working channel for 2.9 mm I=300 mm rigid optics, 12°/30° with Storz connection.

GFRGYN430

Gynko® sheath with 7 FR working channel for 4 mm I=300 mm rigid optics, 12°/30° with Storz connection.

